


# Assista-Rail™ Standard

A standard rail for home beds with accessories to meet changing needs. The Assista-Cane™ is a fixed handle that attaches to the Assista-Rail™ to improve standing and transfer support.

## Features and Benefits


- **Easy Installation** - Rail installs in minutes with no tools. Simply slide the rail between mattress and box spring. Includes safety strap and anti-slide grips to hold rail secure.
- **Sleep Safely** - Effective rail design allows for in-bed mobility & bed queuing (i.e. gives an idea of where the mattress edge is located).
- **Transfer Safely** - By minimizing upper body twisting and reducing wrist strain, the optional Assista-Cane™ provides additional support to reduce the risk of a forward fall.
- **Convenient Storage** - Optional bed pocket provides easy access to personal care items.
- **Dual Rail System** - Install rails on both sides of the bed for added support (conversion kit available).




Assista-Rail™  
Deluxe



Assista-Rail™  
w/ Assista-Cane™



Assista-Rail™  
w/ Bed Pocket



Assista-Rail™  
Standard



Assista-Rail™  
w/ Bed Pocket



TOP MATTRESS REMOVED

## Technical Specifications:

Assista-Rail™	Heavy gauge steel, 7/8" / 2.25 cm diameter. High temperature baked powder coat gloss white enamel. Extends under mattress 24" / 61 cm, up from the mattress 16" / 41 cm and along the side of the bed 18" / 46 cm.
Grip Material	Grey, NPVC non-latex grip. 27" / 68.5 cm long x 1.25" / 3 cm diameter.
Traction Grip Material	Grey, premium quality, non-latex grip. 4.5" / 11 cm long x .84" / 2 cm diameter.
Safety Strap	Strap 1" / 2.5 cm x 100" / 254 cm length with side release buckle.
Assista-Cane™	Heavy gauge steel, 7/8" / 2.25 cm diameter. High temperature baked powder coat gloss white enamel. 31" / 79 cm tall x 7.75" / 20 cm, with grey PVC end cap & rubber cane tip.
Grip Material	Grey, NPVC non-latex grip. 11.5" / 29 cm long x 1.25" / 3 cm diameter.
Bed Pocket	White 420 Denier Nylon, polypropylene strap with velcro attachment.
Weight Capacity	250lbs / 113 kg.
Assistance Level	

**ENTRAPMENT WARNING:** Due to the potential risks associated with head / limb entrapment, the Assista-Rail™ is not recommended for users with the following conditions: erratic body movements, altered mental status, pre-existing conditions such as confusion, restlessness and dementia.